


# PANDA



## How to Bind the Device with Your Account?

## How to Bind the Device with Your Account?

- Step 1: <https://en.freqtek.com:8081> (Enter this website)




**click**

## How to Bind the Device with Your Account?

### Step 2: enter the register slides

#### Register an organization account



The screenshot shows a registration form titled "Register an organization account". At the top, there are two tabs: "clinic" (selected) and "Mechanic / deale". The form contains the following fields:

- \* Name: Please enter full name
- \* phone: +86 Please enter an 11-digit mobile phone number
- \* password: Please enter the password (with a visibility toggle icon)
- \* Confirm: Please enter the password again
- email: Please enter email
- \* Code: Please enter the verification code (with a WVRM logo)

At the bottom, there is a checkbox labeled "I have read [agreement](#)". Below the checkbox are two buttons: "return" and "register".

### Step 3: Prepare these registration information.

(Note: each type of registration cannot use the same phone number and email address, it must be an unregistered phone number and email address)

1. Full name of lab/mechanic/dealer accounts name:
2. Phone number: 15082804444(case)

The phone number requires 11 digits and is in the Chinese number format, (+86...) does not affect the login and use, you can just fill in your phone number and ignore the error note below. Also, the email can also log in and retrieve the password.


Later, it will be possible to register all country phone numbers.

3. email :
4. The password you want to set:
5. Equipment number:ISO....

## How to Bind the Device with Your Account?

### Step 4: Register Devices

Click the “Device” , and then the “Register device” to bind the devices under you account.



The screenshot shows the FREQTY CLOUD interface for 'valerie's clinic'. The left sidebar contains navigation items: Case, Device, Partner, Account, and Dustbin. The 'Device' item is highlighted with a red box and labeled 'Step 1'. The main content area displays a table with columns: number, device status, Device Name, Device ID, Activation time, Expire date, permissions, and Oper. A search bar with 'Device name' and a search icon is present. The 'Register device' button is highlighted with a red box and labeled 'Step 2'. The table currently shows '暂无数据' (No data).


**Note: when the device is bound, it is considered as device activation and the warranty will start with activation.**

## How to Bind the Device with Your Account?

Fill in the information to bind your account.

Note: The device code can only be provided by us according to the Device ID you provided. So before the registration, please give us the device ID and we will tell you the Device Code. The Device ID is on your battery of the device, and you also can check it on the adapter .

That is where you can find the Device ID.



# PANDA



How Can the Clinic  
Send Data to Lab?

## How can the Clinic Send Data to Lab

### Step 1: Cloud platform settings

Register another account for the Lab, and you can use the phone number and email address of the lab to register or register with your own name and email and provide the account to the lab or factory, just to receive data.

(please be noticed when you register an account to received the data, you register under the Mechanic/dealer as the picture shows.)

If you have successfully registered the lab account, you can provide the account to lab or factory to log in and download the data you uploaded after scanning.

Register an organization account

clinic

Mechanic / dealer

\* Name

\* phone

\* password

\* Confirm

email

\* Code


I have read [agreement](#)

return register

## How can the Clinic Send Data to Lab

**Clinic ID:** Log in the clinic account and click **transmission** button of the according device, and add lab information (clinic can send data scanned by this device to lab).

**Lab ID:**no need to click transmission.




The screenshot shows the FREQTY CLOUD interface. The top navigation bar includes the logo, the user name 'ClinicStrayTannlegesenter', a home icon, and a language dropdown set to 'English'. The left sidebar contains navigation options: Case, Device (selected), Partner, Account, and Dustbin. The main content area displays a table of devices with the following columns: number, device status, Device Name, Device ID, Activation time, Expire date, permissions, and Oper. A search bar and a 'Register device' button are located at the top right of the table area. The 'Transmission' button in the 'Oper' column of the first device row is highlighted with a red box.

number	device status	Device Name	Device ID	Activation time	Expire date	permissions	Oper
1	other	PandaP2	IS027E51655ZY	2021-12-13 19:15:54	permanent	user	<a href="#">Add user</a> <a href="#">Transmission</a>



# How can the Clinic Send Data to Lab

Fill in phone number of the lab account.






The screenshot shows the 'FREQTY CLOUD' interface for 'UltradentalindoMandiri'. A 'Data transfer management' dialog box is open, displaying the following information:

- Device ID: IS027E51649ZY
- \* Receiver: Inquire about technical factory (mobile phone n
- Case type: repair, Planting, Orthodontics
- add button
- Table with columns: Case type, Customer name, operation

Case type	Customer name	operation
Planting	Deltalaboratory	delete
repair	Deltalaboratory	delete
Orthodontics	Deltalaboratory	delete

# How can the Clinic Send Data to Lab

Then this lab will show as the partner of this account, and you can upload data to this lab of factory after scanning cases.




# PANDA



How can the Distributor  
Manage the Devices?

## How can the Distributor Manage the Devices?



If you are a distributor, you can bind all devices under your account and then add or authorize users (customers). In this way you can control the user's time of using the scanner.

**Note 1:** When add or authorize a user (customer), you must make sure that the user (customer) also registers a cloud platform account in advance before you can add it successfully.

**Note 2:** When a distributor bind the devices under his account, and add the customers as users, customers cannot see the device list in their platform accounts and cannot set transmission relations neither. They can only see the device list in Appstudio accounts.

**Note 3:** When a distributor bind the devices under his account, and authorize the customers, customers can see the device list and set transmission relations.

## How can the Distributor Manage the Devices?



The screenshot shows the FREQTY CLOUD interface. The top navigation bar includes the logo, the text "s clinic", a language dropdown set to "English", a notification bell, a refresh icon, and a user profile icon labeled "oqit9vtc". The left sidebar contains menu items: "Case", "Device" (highlighted), "Partner", "Account", and "Dustbin". The main content area is titled "Device" and features a search bar with the placeholder "Device Name/Device Number", a search button, a refresh button, and a "Bind device" button. Below this is a table with the following columns: "Number", "Name", "Number", "Activation Time", "Expiration", "permissions", and "operate". A single row is visible with the following data: "1", "PandaP2", a blurred number, "2022-02-14 14:47:29", "permanent", "user", and "operate". Under the "operate" column, there are two buttons: "Add user" (highlighted with a red box) and "Transmission".

Number	Name	Number	Activation Time	Expiration	permissions	operate
1	PandaP2	[blurred]	2022-02-14 14:47:29	permanent	user	<a href="#">Add user</a> <a href="#">Transmission</a>

As a distributor, you can bind all devices under your account, and then click “add user” to authorize the customer using this scanner. You can set the expired time for this scanner.


This is generally recommended when some customers need the payment by installments.

# How can the Distributor Manage the Devices?

### Add user


Device Name PandaP2

\* user

\* **Expire date**  

Remark

Expiration time	Customer name	operation
No Data		



\* Expire date

Remark

Expiration

- permanent
- 7 days
- 15 days
- 1 month
- 3 months
- 6 months
- 1 year
- 5 year


No Data


# PANDA



How to Use the  
Scanner App?

# How to Use the Scanner App




- Step1 : Double-click on the Appstudio  and open it, then log in with the account you registered on our website.  
( Note: Use the account that the devices has been bound, not the lab account. )



# How to Use the Scanner App

**Step 2:** Check if there are devices under “Device” menu. If so then the devices are bound successfully.




The screenshot shows the FREQTY CLOUD interface. On the left, there is a sidebar with three menu items: 'start-up', 'Device', and 'Application'. The 'Device' menu item is highlighted with a red box. The main content area displays a table titled 'Authorization Device' with three columns: 'Device Number', 'Device Name', and 'Authority deadline'. Two rows of data are shown, both highlighted with a red box. The first row has Device Number 'IS027E40511ZY', Device Name 'PandaP2', and Authority deadline 'permanent'. The second row has Device Number 'IS027E52770ZY', Device Name 'PandaP2', and Authority deadline 'permanent'. In the top right corner, there is a 'Refresh the authority' button and a user profile icon labeled 'egzlnsv2'.

Device Number	Device Name	Authority deadline
IS027E40511ZY	PandaP2	permanent
IS027E52770ZY	PandaP2	permanent

## How to Use the Scanner App



**Step 3:** Then double click the Panda P2 App and open it. If it shows the page like this, then we have successfully installed all the software and you can use it right now!

A screenshot of the PandaP2 application interface. The top navigation bar includes "FREQTY" and several menu items: "Orders", "Planning", "Acquisition", "Processor", and "Storage". A left sidebar contains icons for "Create", "Admin", "Import(O)", "Practice(E)", and "Recent". The main content area displays a "Case Information" form with fields for "Patient Name" (a dropdown menu), "ID" (containing "220125001"), "Dentist" (a dropdown menu), and "Clinic" (a text input field). At the bottom of the form are "Exit" and "Next" buttons.