





PANDA SCANNER **ALL OVER THE WORLD**

www.panda-scanner.com info@panda-scanner.com

Chengdu | Ziyang | Suzhou | Ningbo | Hong Kong | Brazil | USA



Manda *
Scanner







Smaller and Lighter

Only 138g 221x27x25mm



Faster and More Accurate

l min ultra-fast scan 7 μm ultra-high accuracy



Plug and Play

No power adapter needed



New Designed Tip

White Silver 50-100 autoclaving cycles
More ergonomic



New Collection

The art of intraoral scanners

- Crystal White
- Radiant Silver
- Sparkly Purple
- Ice Blue
- Lightsaber Red
- Champagne Gold
- Laguna Blue





Ergonomic

Efficient

Innovative




SOFTWARE FEATURES



PREDITY © B. Order S Parring X Appetition O Processor O Borrige O Company Description of the Company of the Co







Undercut Check

Auxiliary confirmation in place, colour gradient indicates undercut, accurate prediction, and improves the communication efficiency between doctors and technicians.

Margin Line Drawing

Draw accurate and clear margin lines, linking doctors and technicians to reduce rework.

Pre-preparation Scan

Obtain pre-preparation tooth state, retain bite data, and reshape the original good teeth.

Bite Check


Color gradient indicates distance, accurately judges the restoration distance and assists in improving the effect of tooth preparation.


Lock Data

Locks the tooth position to prevent data loss and reject excess data interference.











Color Switch

True color jaws, clearly restore details, one click to switch colors, and check model data comprehensively.

Al Intelligent Scan

Easily identify soft and hard tissues, eliminate redundant data, and make accurate dental models.

Orthodontic Simulation

Generate simulation animations to help doctors and patients communicate, share quickly, and promote orders.


Disease Course Monitoring

Full-cycle diagnosis and treatment monitoring, comparative observation of model differences, simultaneous analysis, comprehensive upgrade of diagnosis and research.

Cuff Tool

Clear cuff data, avoid tissue interference and lock automatically, significantly improving data reliability.

APPLICATIONS



SPECIFICATIONS

PRODUCT INFORMATION

Category	PANDA smart
Anti Fog System	Intelligent Heating
Accuracy	≤7µm
Precision	≤5µm
Handpiece Weight	138g (192g with cable)
Handpiece Dimension	221x27x25mm
Handpiece Enclosure	High-strength Aviation Aluminum Alloy
Scan Field	18x16mm
Scan Depth	0-20mm
True Color	Full HD
Tips	x6 (Autoclavable)
Tip Autoclave Cycle	50-100 Times
Data Interface	USB3.0
Cable Length (m)	1.85m



SOFTWARE

Category	PANDA smart
Types of Scanning	Orthodontics, Implant, Restoration
Image Capture	Yes
3D Data Preview	Yes
Output Files	STL, PLY, PTY (Encrypted for user)
Cloud Service	Yes
Cloud Service Price	5GB/Free
Direct Compatibility with Dental CAD	Exocad*
Main Image Processing Tools	Trim/Edit/Recover/Lock, etc.
Margin Line Drawing	Yes
Occlusal Match	Yes
Undercut	Yes
Data Comparison	Yes
Al Technology	Auto Tissue Remove, Overlap Correct, etc.
Software License	Free
Software Upgrade	Free
Transfer Data with API	Available

PC REQUIREMENTS

Category	PANDA smart
CPU	Intel i7-12700H/Intel i7-11800H or above
RAM	16G/32G
Hard Disk	SSD 512G or above
GPU	RTX2060/RTX3060 or above
Operating System	Windows10/11 64bit
Display	Resolution 1920×1080 and above

^{*}Check the PC Recommendation Catalog on our website for more detail.

WARRANTY

Category	PANDA smart
Standard Warranty	2 years
Extended Warranty	2 years more available

^{*}All names, logos, brands and other trademarks featured or referred to within this brochure are the property of their respective trademark holders.

CUSTOMER SUPPORT

Email: support@freqty.com
Toll Free Number (USA): 888-855-2562
WhatsApp (Global): +86 18108286761

